

Manual-Active8

Super Concentrated Manual Enzymatic Detergent for Cleaning Reusable Medical Instruments & Devices

Product Properties:

- Enzymatic Detergent.
- Great cleaning power.
- Economical.
- Removes any organic traces e.g., blood, protein, fat, saliva, chemical residues.

Material Compatibility:

Manual-Active8 is designed for the manual cleaning of flexible endoscopes and reusable surgical instruments.

*Stainless Steel *Chromium Plated Material *Synthetic Materials *Rubber
*Latex *Glass *Ceramics *Copper *Brass *Metal Alloys

When used correctly in the correct dilutions, no material damage is to be expected.



Areas of Application:

- For use with manual cleaning of reusable surgical instrumentation.

Composition:

Manual-Active8 contains amongst other ingredients:

Monopropylene Glycol CAS No. 57-55-6, Propan 2-OL CAS No 67-63-0, Subtilisin CAS No 9014-01-1
pH 7.5-8.7 typical

Directions of Use:

Use to remove residues and deposits found on and in flexible devices. **Dilution Rate:** Inject 4ml per litre of wash water. **Temperature Range:** Between 27-35°C. **Minimum Contact Time:** 2 Minutes. **Shelf Life:** 24 Months.

Ordering Details:

Product Code:	1039	10591	1059	1059U	105910
Pack Size:	6x500ml	12x1L	2x5L	UN 2x5L	10L