

## Delta E.S

### Super Concentrated Enzymatic Alkaline Instrument Washing Compound

#### Product Properties:

- Mildly Alkaline.
- Phosphate free non-foaming washing compound.
- Designed for the cleaning of re-usable instrumentation in the automated phase.
- Economical.
- Suitable for all makes of automated washers.
- Excellent cleaner in all water types.
- Excellent material compatibility.
- Will not cause corrosion or oxidation.
- Aqueous solution dispensed in the usual way.



#### Material Compatibility:

Delta E.S is designed for the automated cleaning of surgical instruments.

\*Stainless Steel \*Chromium Plated Materials \*Synthetic Materials \*Rubber  
\*Latex \*Glass \*Ceramics \*Copper \*Brass \*Metal Alloys

When used correctly in the correct dilutions, no material damage is to be expected.

#### Areas of Application:

- For use in washer disinfectors.

#### Composition:

Delta E.S contains amongst other ingredients:

Monopropylene Glycol CAS No 57-55-6, Tetrasodium Bisphosphonate CAS No 3794-83-0,

Subtilisin CAS No 9014-01-1

pH 10-11 typical

#### Directions of Use:

DO NOT USE BY HAND. When changing containers, wear rubber gloves, goggles, and plastic apron. Place feed tube into bottom of container and check automatic feed equipment is switched on.

**Dilution Rate:** Inject 2ml per litre wash water at a maximum working temperature of 65°C, or 4ml per litre wash water at a maximum working temperature of 60°C. **Temperature Range:** 60°C-65°C. **Minimum Contact Time:** 2 Minutes. **Shelf Life:** 24 Months.

#### Ordering Details:

Product Code:	107021	107021U	107020	107022
Pack Size:	2x5L	UN 2x5L	10L	20L