

P.D.A. 51

Alkaline Hospital Instrument Cleaner Concentrate for use in Washer Disinfectors, Ultrasonic Cleaners, Sonic Irrigators & Manual Cleaning

Product Properties:

- Concentrated alkaline detergent.
- Designed for automated or manual cleaning of reusable medical instruments and devices, such as surgical instruments.
- Safe for use on soft metals such as aluminium and anodised aluminium.
- Economical.

Material Compatibility:

P.D.A. 51 is designed for the automated or manual cleaning of hospital instruments.

*Stainless Steel *Aluminium *Anodised Aluminium *Chromium Plated Materials *Synthetic Materials *Rubber *Latex *Glass *Ceramics *Copper *Brass *Metal Alloys

When used correctly in the correct dilutions, no material damage is to be expected.



Areas of Application:

- For use in the cleaning of reusable medical instruments and devices.

Composition:

P.D.A. 51 contains amongst other ingredients:

Potassium polyphosphate aqueous solution, Disodium Metasilicate CAS No 6834-92-0,

Subtilisin CAS No 9014-01-1

pH 13-14 typical

Directions of Use:

When changing containers, wear rubber gloves, goggles and plastic apron. Place feed tube into bottom of container and check automatic feed equipment is switched on.

Dilution Rate: Manual Cleaning: Use 3-6ml per litre of tepid water.

Automated Cleaning: Inject 3-6ml per litre of water at a working temperature of 15-70°C.

Temperature Range: Between 15-70°C. **Shelf Life:** 24 months.

****Do NOT add any other chemicals such as bleach or detergents, as these may destroy the effectiveness of this product****

****Consult the manufacturers recommended reprocessing instructions to determine cleaning agent compatibility and any additional procedures necessary to clean their medical device****

Ordering Details:

Product Code:	PDA51/5	PDA51/5U	PDA51/10	PDA51/20
Pack Size:	2x5L	UN 2x5L	10L	20L