

Process pH Plus

Alkaline Instrument Cleaner Concentrate

Product Properties:

- Powerful cleaner designed for hand washing of instruments.
- Great cleaning power.
- Economical.
- Removes any organic traces e.g., blood, protein, fat, saliva, chemical residues.
- Complimented with a very pleasant and fresh smelling fragrance.

Material Compatibility:

Process pH Plus is designed for the manual cleaning of reusable surgical instruments.

*Stainless Steel *Chromium Plated Material *Synthetic Materials *Rubber
*Latex *Glass *Ceramics *Copper *Brass *Metal Alloys

When used correctly in the correct dilutions, no material damage is to be expected.



Areas of Application:

- For use with manual cleaning of reusable surgical instrumentation.
- NOT TO BE USED FOR ENDOSCOPE REPROCESSING

Composition:

Process pH Plus contains amongst other ingredients:
Disodium Cocoiminodipropionate, Disodium Metasilicate CAS No 6834-92-0
pH 12-14 typical

Directions of Use:

Use for the hand washing of instruments for removal of bioburden and body fat deposits.

Dilution Rate: For general cleaning: Dilute 3ml to 5 litres of wash water.

For heavily soiled items: Dilute 6ml to 5 litres of wash water.

Temperature Range: 27°C-45°C. **Minimum Contact Time:** 2 Minutes. **Shelf Life:** 24 months.

Ordering Details:

Product Code:	3063TWIN	3064	3063	3063U	3068	3069
Pack Size:	6x500ml	12x1L	2x5L	UN 2x5L	10L	20L