

Progenica

Ultra Concentrated Alkaline Washing Compound for use in Washer Disinfectors & Ultrasonic Machines

Product Properties:

- Mildly Alkaline.
- Phosphate free non-foaming washing compound.
- Designed for the cleaning of reusable instrumentation in the automated phase.
- Economical.
- Suitable for all makes of automated washers.
- Excellent cleaner in all water types.
- Aqueous solution dispensed in the usual way.

Material Compatibility:

Progenica is designed for the automated cleaning of surgical instruments.

*Stainless Steel *Chromium Plated Materials *Synthetic Materials *Rubber
*Latex *Glass *Ceramics *Copper *Brass *Metal Alloys

When used correctly in the correct dilutions, no material damage is to be expected.



Areas of Application:

- For use in washer disinfectors.
- Can also be used in ultrasonic machines.

Composition:

Progenica contains amongst other ingredients:

Hydroxyethylidenediphosphonic Acid CAS No 29329-71-3, Disodium Metasilicate CAS No 6834-92-0, Sodium Hydroxide CAS No 1310-73-2

pH 12-14 typical

Directions of Use:

DO NOT USE BY HAND. When changing containers, wear rubber gloves, goggles and plastic apron. Place feed tube into bottom of container and check automatic feed equipment is switched on.

Dilution Rate: Inject 2ml per litre wash water at a maximum working temperature of 65°C or 4ml per litre at a maximum working temperature of 60°C. **Temperature Range:** Maximum working temperature of 65°C.

Minimum Contact Time: 2 Minutes. **Shelf Life:** 18 Months.

Ordering Details:

Product Code:	1098	1095	1095U	1096	1097	2095	5095	101095
Pack Size:	12x1L	2x5L	UN 2x5L	10L	20L	200L	500L	1000L