

Sequel

Acid Free Rinse Agent for Instrument Washing Machines

Product Properties:

- Acid free rinse aid.
- Will aid drying process of reusable surgical instrumentation.
- Low Foaming.

Material Compatibility:

Sequel is designed for the automated drying of surgical instrumentation.

*Stainless Steel *Aluminium *Anodised Aluminium *Chromium Plated
Materials *Synthetic Materials *Rubber *Latex *Glass *Ceramics *Copper
*Brass *Metal Alloys

When used correctly in the correct dilutions no material damage is to be expected.



Areas of Application:

- For use in the final rinse phase during the decontamination process.
- For use in automated dispensing systems.

Composition:

Sequel contains amongst other ingredients:
Isopropyl Alcohol CAS No. 67-63-0
pH 7-8 typical

Directions of Use:

Sequel is for use through automatic dispensing systems only.

Dilution Rate: Inject 1-2ml per 5 litres of final rinse water.

Temperature Range: Between 82-95°C.

Ordering Details:

Product Code:	10741	1074	1075	1076	1085	2074	2075	101074
Pack Size:	12x1L	2x5L	UN 2x5L	10L	20L	200L	500L	1000L