

Ultra-Active8

Machine Enzymatic Detergent for Cleaning Reusable Medical Instruments & Devices

Product Properties:

- Enzymatic detergent.
- For use in ultrasonic washers.
- Great cleaning power.
- Economical.
- Removes any organic traces e.g., blood, protein, fat, saliva, chemical residues.

Material Compatibility:

Ultra-Active8 is designed for the automated cleaning of surgical instruments.

*Stainless Steel *Chromium Plated Material *Synthetic Materials *Rubber
*Latex *Glass *Ceramics *Copper *Brass *Metal Alloys

When used correctly in the correct dilutions, no material damage is to be expected.



Areas of Application:

- For use in ultrasonic washers.
- For use in automated washer disinfectors.

Composition:

Ultra-Active8 contains amongst other ingredients:

Monopropylene Glycol CAS No 57-55-6, Propan 2-OL CAS No 67-63-0, Subtilisin CAS No 9014-01-1
pH 7.5-8.7 typical

Directions of Use:

When changing containers, wear rubber gloves, goggles, and plastic apron. Place feed tube into bottom of container and check automatic feed equipment is switched on. **Dilution Rate:** Inject 4ml per litre of wash water. **Temperature Range:** Between 27-60°C. **Minimum Contact Time:** 2 Minutes. **Shelf Life:** 24 Months.

Ordering Details:

Product Code:	10510	1051	1051U	1053	1053200
Pack Size:	6x500ml	2x5L	UN 2x5L	10L	200L