

Flexi-Clenz Bedside Kit

Non-Enzymatic Bedside Response Kit

Product Properties:

- Versatile in application.
- Economical.
- Removes any organic traces e.g., blood, protein, fat, saliva, chemical residues found on and in flexible and rigid devices.
- Complies with BSG guidelines.
- Suitable for all makes and manufacturer of flexible and rigid endoscopes.
- Lightly foaming.



Material Compatibility:

Flexi-Clenz Bedside Kit is designed for the manual cleaning of flexible and rigid endoscopes and reusable surgical instruments.

*Stainless Steel *Aluminium *Anodised Aluminium *Chromium Plated Material *Synthetic Materials
*Rubber *Latex *Glass *Ceramics *Copper *Brass *Metal Alloys

Areas of Application:

- For use with flexible and rigid endoscopes.

Composition:

Flexi-Clenz Bedside Kit contains amongst other ingredients:
Propan-2-OL CAS 67-63-0, Glycerine (Glycerol) CAS No 56-81-5
pH 6-7.5 typical

Directions of Use:

Use to remove residues and deposits found on and in flexible and rigid endoscopes.

Dilution Rate: Ready to Use.

Ordering Details:

Product Code:	1031	1033
Pack Size:	27x250ml	27x500ml